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Universal submicron hybrid composite that achieves high aesthetic results using only a small range of shades.

•• INDICATIONS

- Highly esthetic restorations of I, II, III, IV and V cavities;
- Onlay and Inlay;
- “Sandwich” technique;
- Fiber Splinting;
- Crown Build Up.

•• INSTRUCTION

Cavities are prepared in the conventional manner. In very deep or large restorations and in cases of near pulp exposure the use of a calcium hydroxide liner is recommended e.g. Jen-Line LCF.

•• RESTORATING OF ANTERIOR GROUPS OF TEETH

The use of a mylar strip is recommended for interproximal separation and as a matrix. Etch the enamel surrounding the cavity including undercuts with enamel conditioner [Phospho-Jen AS Gel] for 15 – 20 sec., wash and dry. Apply Bonding Agent [Jen-Unibond is recommended] over the enamel margins and cure for 10 seconds. Place restorative material with a tapping movement and cure. In very deep cavities more than one layer may be required. Shape with appropriate finishing diamond or multi-fluted carbide burs, contour proximal surfaces with finishing strips. Check occlusion with a thin articulating paper, carefully adjust occlusion by removing material with a fine polishing diamond or stone. Finish with discs and rubber cups.

•• RESTORATING OF POSTERIOR GROUPS OF TEETH

Place a mylar strip or matrix band for interproximal separation and fix it with wedges. Adapt the band to seal the gingival area to avoid overhangs. Note: the matrix may be placed following the enamel etching if preferred. Apply Bonding Agent [Jen-Unibond is recommended] over the enamel margins and cure for 10 seconds. To aid in adaptation, the first 1,5 mm layer of Jen-Radiance may be placed and adapted to the proximal box. A condensing instrument can be used to adapt the material to all of the internal cavity aspects and cure it. The next portions of Jen-Radiance may be placed in layers less than 2,5 mm. Cure each increment separately. Shape with appropriate finishing diamond or multi-fluted carbide burs, contour proximal surfaces with finishing strips. **Attention!** Jen-Radiance Opaque shades should be applied in layers no more than 1 mm thick.

•• CAUTIONS

Jen-Radiance contains methacrylic resins that may cause an allergic reaction by skin contact in certain individuals. Avoid long or repeating contact of not-polymer

•• PRECAUTIONARY MEASURES

1. During the work with the material use special glasses, gloves, clothes and mask Goggles are recommended for the patients also.

2. Observe the security measures provided during the work with powdery, siliceous substances. Use goggles, a mask and gloves.

•• CONTRAINDICATIONS

Allergic reactions in the anamnesis of the patient especially on methacrylic resins or on any other component of dental materials.

•• COLLATERAL REACTIONS

The product can cause irritation of eyes, skin and mucous membranes. (see. CAUTION section).

•• INTERACTION WITH OTHER DENTAL MATERIALS

Use of eugenol-containing materials in a combination with Jen-Radiance is contraindicated. Evgenol-containing dental materials can have negative impact on ability to an polymerization.

•• WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty. The customer is responsible for determination of suitability of this product for user’s application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

•• LIMITATION OF LIABILITY

The manufacturer’s liability is limited by only cases stipulated by direct legislation of the country.

•• STORAGE

· Inappropriate storage conditions will reduce shelf life and can lead to deterioration of properties of the material. Don’t expose this product to direct sunlight. Store a material in a dry place.

· Store this product at the temperature of 4 – 25 °C [39,2 – 77 °F].

· DO NOT FREEZE!

· Shelf life of Jen-Radiance is 3 years.

•• RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

•• PACKING

Jen-Radiance is packed in syringes on 4 g.

· Translucent (incisal) shades: I (Incisal - Enamel Normal), IOP (Incisal Opaline - Enamel White Opalescent);

· Universal dentin shades: A1, A2, A3, A3.5, B1, B2, B3, C2, D2, BW (Bleach White), XBW (Extra Bleach White);

· Opaque shades: A1-0, A2-0, A3-0, A3.5-0, UO (Universal Opaquer), SWOP (Super White Opaquer).

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